

**⚠ IMPORTANT: REPORT ANY SHIPPING DAMAGE IMMEDIATELY:**  
**Inform shipper of any damages - Leave carton intact. Leave the equipment in the receiving area until inspection is complete.**

**⚠ CAUTION: POSSIBLE EQUIPMENT DAMAGE:**  
**This crate contains fragile expensive medical equipment. Uncrate and handle carefully. If after uncrating the equipment you find any damage (no matter how slight), show it to your supervisor.**

**⚠ WARNING: PERSONAL INJURY HAZZARD:**  
**When cutting bands, always use a tool specifically designed for that purpose. This will help avoid personal injuries frequently incurred when bands are cut and tension is released.**

**IMPORTANT:** Follow each step in the order shown in these instructions.

### **UNPACKING INSTRUCTIONS:**

Your Hausted equipment has been carefully packed at our manufacturing plant to ensure safe shipment to your medical facility. There are several procedures you must follow to put your new equipment in service. These procedures only take a few minutes to complete and are required to ensure proper operation of the equipment.

1. Cut the two bands around the shipping carton.
2. Remove the top half of the carton and cut one side of the bottom half.
3. Ensure the fifth wheel / power drive wheel clears the edge of the skid when removing product (where applicable), then remove the equipment from the carton.
4. Check to see if all features of the equipment work properly. If all the features work, advance to step 5. If any of the features of the equipment do not work properly, call GF Health Products, Inc. for service at: 770.368.4700.
5. Clean the equipment using mild detergent to remove any dirt accumulated during shipment, and place the equipment into service.