

# FLUORO-TRACK STRETCHER



FLUOROSCOPY AND C-ARM ACCESS STRETCHER

## 493RPCST FLUORO-TRACK STRETCHER

Designed to offer the security of a treatment table with the mobility of a stretcher. The Fluoro-Track Stretcher is designed for use with C- Arms for fluoroscopy-guided procedures. Its unique design, with a proprietary sliding top, enables a broader range of imaging capacity without repositioning the patient.

### SLIDING TOP:

- Easy to operate sliding top enables a broad range of imaging capacity without repositioning the patient, enabling shorter procedure times
- Full length, radiolucent tops translate to either end providing excellent patient access
- 18 inch total range top slides 12 inches (30.5 cm) from center position toward foot end and 6 inches (15.2 cm) toward head end




### FEATURES:

- Low profile trimmed base creates an excellent window for C-arm access from both sides
- Designed to be a procedure table with the mobility of a stretcher, it provides improved patient load ratio for easier transport reducing risk of patient and staff injury with patient transfers
- Dual side foot pedal hydraulic controls for easy height adjustment
- Fold away side arm rails with built-in push-steer handles
- Dual column base provides superior support for the patient
- 8 inch Tente® casters facilitate maneuverability
- 3 inch comfort pads

★ See Chart for Options and Specifications

493 Product Info



### Specifications: 493RPCST

	inch	cm†
Overall Width:	29.75	76
Overall Length:	80.25	204
Overall Height: range - minimum	23	58.4
Overall Height: range - maximum	36	91
patient surface Width:	24	61
patient surface Length:	76	193
Retracto™ Rail Length: upright position	64	162.5
Retracto™ Rail Height: upright position	15	38
Retracto™ Rail Height: mid-position	12.75	32.4
C-Arm Clearance	23	58.5
Radiolucent Pad: full length	3	7.6
Tente® Casters: easy rolling	8	20.3
	degree of movement™	
Backrest: movement	0 - 80	
Trendelenburg: standard	0 - 15	
Trendelenburg: reverse	0 - 15	
Weight Capacity: EVENLY DISTRIBUTED	lb	kg†
	500	227
Shipping Weight:	364	165

### Features:

Operation: Hydraulic

■ Standard ■ Optional — Does Not Apply  
† standard rounding method used to show metric values as whole numbers

Product Hyperlinks Enabled



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2014, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.