



INTRODUCING THE HSS SERIES SURGICAL STOOL MEDICAL SEATING



HSS500-W WEDGE
SEAT STYLE

HSS500-S SADDLE
SEAT STYLE
SHOWN RIGHT



HSS500-WD WATER DROP
SEAT STYLE

HSS500-RWD REVERSE
WATER DROP
SEAT STYLE
ALSO AVAILABE



BUILT FOR
PREFORMANCE
AND COMFORT

HAUSTED SURGICAL STOOL

IDEAL FOR EYE, ENT AND PLASTIC PROCEDURES



- Lightweight, maneuverable unit fits easily through doorways
- Three-Spindle Base and four dual Tente® casters for greater stability
- Dual-density foam seat with Dura-Gel® insert for enhanced comfort
- 360° Arm Rest Rotation, 180° arm rest support rotation
- Foot-controlled hydraulics

Adjustable Contoured Backrest

[HSS500 Product Info](#)



360° Rotating Urethane Arm Rest for Comfortable Arm Positioning

Dual-density Foam Seat with **DURA-GEL®** Insert for Enhanced Comfort



180° Arm Rest Support for more Stable Arm Positioning

Adjustable Height Range
21 - 28 inches (53 - 71 cm)

ABS Thermoformed Hood for Easy Cleaning

OPTIONAL ACCESSORY:
Riser Kit increases Height Range an additional 2 inches (5 cm)
Item No: HSS500RISER



FOOT-CONTROLLED HYDRAULICS with Three-Spindle Base and Four 3 inch Dual Tente® Casters for Greater Stability (8 cm)



Front Locking System

Options:	Seat Width:		Seat Depth:	
	inch	cm	inch	cm
HSS500-S saddle seat	17.00	43.18	17.65	44.83
HSS500-W wedge seat	13.50	34.29	18.00	45.72
HSS500-WD water drop seat	15.75	40.00	18.00	45.72
HSS500-RWD reverse water drop seat	15.75	40.00	18.00	45.72



Specifications: **HSS500-xx**

	inch	cm [†]
Overall Width: [§]	25.50	64.77
Overall Depth: [§]	32.25	81.92
Overall Height: minimum	21.00	53.34
Overall Height: maximum	28.00	180.34
Base Width: [§]	24.00	60.96
Base Depth: [§]	28.50	72.39
Seat Height: minimum	21.50	54.61
Seat Height: maximum	28.50	72.39
Seat Thickness: dual-density foam	3.00	10.16
Backrest Width:	15.56	39.52
Backrest Length:	8.25	20.96
Backrest Height Range: minimum	10.50	26.67
Backrest Height Range: maximum	12.50	31.75
Armrest Width: [§]	4.75	12.07
Armrest Length: [§]	12.00	30.48
Armrest Height: minimum	7.00	17.78
Armrest Height: maximum	15.00	38.10
Caster Diameter:	3.00	7.62

300 POUND¹
WEIGHT CAPACITY
EVENLY DISTRIBUTED

	features	
Caster Type:	Dual Tente®	
	degrees [°]	
Armrest Rotation: movement	360	
Armrest Support: movement	180	
Backrest Tilt: from vertical	20	
	lb	kg [†]
Weight Capacity: ¹ EVENLY DISTRIBUTED	300	136.08
Shipping Weight: [§]	74	33.57

¹ maximum weight capacity - evenly distributed - refer to user manual
[†] standard rounding method used ² variance ±0.5 inch
[§] approximate measurement

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company
+1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH