

HORIZON[®] RETRACTO[®] RAIL HYDRAULIC STRETCHER



DESIGNED FOR PROCEDURE, RECOVERY & TRANSPORT

462 HORIZON® SERIES RETRACTO® RAIL HYDRAULIC STRETCHER

The Retracto™ Rail Hydraulic Stretcher is designed to provide longlasting and consistent performance. These easy-to-clean units support patients during transport in any hospital or surgical environment.

VERSATILITY:

- Streamline patient handling allowing this unit to be used from preinduction to transport to recovery to help improve case turnaround
- Easy to use hydraulic foot pedals to adjust height and Trendelenburg on either side of the stretcher
- Intermediate side rail position offers maximum patient access without compromising patient or staff safety
- Dual accessory wells at all four corners for versatility in adapting accessories
- Full perimeter heavy duty non-marring bumper system
- Seamless thermoformed base cover with storage
- Built-in oxygen tank holder

MOBILITY:

- Three modes of operation: steer, neutral and brake
- Brake and steer levers are accessible at all four corners of the stretcher base
- 8 inch Tente® Casters



462 Product Info



Specifications: 462 RETRACTO

	inch	cm†
Overall Length:	83	210
Overall Height: <i>range - minimum</i>	22.25	56
Overall Height: <i>range - maximum</i>	34.75	88
side rail Length:	60	152
side rail Height:	15	38
Tente® <i>easy rolling</i> Casters:	8	20
	degree of movement*‡	
Cable Free Fowler: <i>pneumatic-assist</i>	0 - 90	
Trendelenburg: <i>standard</i>	0 - 16	
Trendelenburg: <i>reverse</i>	0 - 16	
	lb	kg†
Shipping Weight:	330	150

Features:

Rail Type:	Retracto™
IV Wells:	4 - 1/corner
Accessory Wells:	4 - 1/corner
Brake/Steer Pedals:	4 - 1/corner
Hydraulic Controls:	Foot Pedal - 1/side

426HMCST PICTURED
cover: with arm board accessory
top: in trendelenburg position
right: with back tilted up



	462EMCST		462FMCST		462GMCST		462HMCST	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Overall Width:	30	76	34.5	88	30	76	34.5	88
patient surface Width:	25	64	30	76	25	64	30	76
patient surface Length:	76	193	76	193	76	193	76	193
	feature		feature		feature		feature	
Knee Gatch:	No		No		Yes		Yes	
	lb	kg†	lb	kg†	lb	kg†	lb	kg†
Weight Capacity: <small>EVENLY DISTRIBUTED</small>	500	227	800	363	500	227	800	363

■ Standard ■ Optional — Does Not Apply
 † same specs for both AC110 volt and the DC220 volt versions
 ‡ standard rounding method used to show metric values as whole numbers



GF Health Products, Inc. is an
 Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
 Hausted is a registered trademark of GF Health Products, Inc.

GF1400057/RevF25