



SUPPORT SURFACES



HAUSTED COMPETITOR REPLACEMENT STRETCHER PADS

HSP STANDARD SERIES STRETCHER PADS

The Hausted® Standard Stretcher Pads are for use with popular Stretcher models from many manufacturers. The pads offer exceptional cushioning, support, and pressure redistribution. The fabric is durable and offers low friction for easy patient transfers.

Features:

Standard foam pad height

- 2 layer foam pad for comfort and durability
- Hook attachment on base
- Easy to clean and maintain

Cover Features:

- Herculite® Sure-Chek™ cover material is comfortable and durable
- Infection Control is aided by Antimicrobial properties; Bacteria resistant, anti-fungal, and non-allergenic Flame Resistant
- Fluid Resistant
- Anti-static
- Stain Resistant



600C COMFORT CARE SERIES STRETCHER PADS

The Hausted® Clinical Comfort Stretcher Pads are the best choice when selecting a non-powered support surface. These pads are for use with popular Stretcher models from many manufacturers.

Features

5 inch thick foam pad

- Tri-laminate foam pad design for better comfort and support
- Contour cut grid allows for better patient immersion, and greater airflow
- Better pressure redistribution properties than standard pad offerings
- Hook attachment on base
- Easy to clean and maintain

Cover Features:




- Herculite® Sure-Chek™ Fusion III cover material is a 4-way stretch
- Infection Control is aided by Antimicrobial properties; Bacteria resistant, anti-fungal, and non-allergenic
- Flame Resistant
- Fluid Resistant
- Anti-static
- Stain Resistant



STANDARD SERIES
2-layer foam construction for comfort



600C — COMFORT CARE SERIES
3-layer construction with contour-cut grid pattern provides greater airflow and improves patient immersion



Specifications:


	Standard	HSP-600C-	1007-01	1007-02	1007-03	1007-04	1007-05	1007-06	1007-07	1007-08	1007-09	1007-10
	Clinical Comfort	600C-	100701	100702	100703	100704	100705	100706	100707	100708	100709	100710
		key	inch	cm [†]	inch	cm [†]	inch	cm [†]	inch	cm [†]	inch	cm [†]
Overall Width:	A		76	193.04	76	193.04	76	193.04	76	193.04	75	190.50
Overall Length:	B		26	66.04	29	73.66	31	78.74	30	76.20	27	68.58
Overall Height:	C		4	10.16	4	10.16	4	10.16	4	10.16	4	10.16
Overall Height:	C		5	12.70	5	12.70	5	12.70	5	12.70	5	12.70
Head Taper:	D		9	22.86	9	22.86	9	22.86	7	17.78	8.5	21.59
Foot Taper:	E		4	10.16	4	10.16	4	10.16	4	10.16	8.5	21.59
Foot Taper:	F		4	10.16	4	10.16	4	10.16	4	10.16	4.5	11.43
distance from Foot:	G		26	66.04	26	66.04	26	66.04	26	66.04	10	25.40
distance from Head:	H		—	—	—	—	—	—	—	—	29	73.66
Hook Width [‡] :	J		2	5.08	2	5.08	2	5.08	2	5.08	6	15.24
Hook Length [‡] :	K		24	60.96	24	60.96	24	60.96	24	60.96	4	10.16
Number of Hooks [‡] :			1		1		1		1		2	
hook/loop attachment												
Warranty:			12		12		12		12		12	
in months												
Weight Capacity:			500	226.80	500	226.80	660	299.37	500	226.80	500	226.80
EVENLY DISTRIBUTED			lb	kg [†]	lb	kg [†]	lb	kg [†]	lb	kg [†]	lb	kg [†]

[‡] Location of hook/loop on bottom of pad may not match location of attachment point on stretcher. Make sure pad is securely attached to stretcher to avoid injury.
[†] standard rounding method used to show metric values as whole numbers

STRETCHER FABRICS

The Hausted stretcher pads use fabrics designed to meet the demands in a healthcare environment. Pads are fluid-resistant and resist staining from most contaminants such as medications, oils, and more. The fabric contains antimicrobial additives designed to help protect the fabric against odor-forming bacterial growth and reduces undesirable organic odors.

STANDARD STRETCHER PAD FABRICS



036 BLACK

SURE-CHEK® COMFORT FABRIC

Sure-Chek Comfort is a fluid and moisture-resistant cover known for its durability and easy cleaning. It is a Vinyl laminate material that conforms to the body and allows the inner construction of the stretcher pad to provide comfort and support.

UPGRADED STRETCHER PAD FABRICS



043 SPRUCE

SURE-CHEK® FUSION III®

The **Sure-Chek Fusion III** is engineered for maximum performance. The 4-way stretch, fluid-resistant, urethane cover reduces friction and skin shear while helping to keep the patient cool and comfortable. The fabric provides maximum resistance against spills, staining, and odors. The fabric incorporates antimicrobial additives designed to help protect the fabric against odor-forming bacterial growth to reduce undesirable organic odors.

Fits These Models:

- Advantage: 1001/1009
- Advantage: 1501/1509
- Atlas: 660
- M series: — 26" Litter
- M series: — 30" Litter
- Renaissance: 1210/1231
- Renaissance: 1710/1731
- Transport: 720/721 — 24" Litter
- Transport: 720/721 — 26" Litter
- Transport: 737
- Trauma: 1002/1020
- Trauma: 1027

Fits These Models:

- Electric: (P8020)
- GPS: 880/881/882/883
- GPS: 885/886
- Procedural: (P8000)
- Surgical: (P8010)
- Transport: (P8005) — 26" Patient Surface
- Trauma: (P8040)

Our Part N°:	Stryker®											
xxxx-100701	■			■		■				■	■	■
xxxx-100702		■					■					
xxxx-100703			■									
xxxx-100704				■								
xxxx-100705							■					
xxxx-100706										■		

Our Part N°:	Hill-Rom®						
xxxx-100707	■						
xxxx-100708		■					
xxxx-100708			■				
xxxx-100709				■		■	■
xxxx-100710					■		
xxxx-100706							

■ Our Part N° fits Stryker/Hill-Rom Model N° Shown

— to order replace the 'xxxx' in the our part number with 'HSP' (standard pad) or '600C' (comfort pad)

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.hausted.com | www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved. Hausted is a registered trademark of GF Health Products, Inc. Stryker is a registered trademark of Stryker Corporation. Hill-Rom is a registered trademark of Hill-Rom Holdings, Inc.



GF160003 Rev D25