

# REPLACEMENT TOPS



## REPLACEMENT TOPS

Make your exam tables look new at a fraction of the cost. All components within in the replacement top are all new including wood, foam, vinyl and hardware.

- Industry standard colors *plus*
- Replacement Tops fit leading brand tables



Specifications:	100   104		204 <sup>1</sup>   223 <sup>2</sup>		304		404   405		420   460	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
Overall Width:	27.50	69.85	27.50	69.85	27.50	69.85	27.50	69.85	27.50	69.85
Overall Length:	55.00	139.70	55.00	139.70	55.00	139.70	60.00	152.40	60.00	152.40
Footrest Width:	20.00	50.80	19.00	48.26	20.00	50.80	20.50	52.07	19.00	48.26
Footrest Length:	18.00	45.72	18.00	45.72	18.00	45.72	17.50	44.45	18.00	45.72
Shipping Weight:	55	24.95	55	24.95	55	24.95	35	15.88	55	24.95

Features:	100   104	204 <sup>1</sup>   223 <sup>2</sup>	304	404   405	420   460
Velcro® Attachment:	—	—	—	included	—
Paper Roll Holder:	included	—	—	—	included
Mounting Hardware:	included	—	—	—	included
Lift Handles:	—	not included <sup>3</sup>	—	—	—
Top Replacement Knobs:	—	not included <sup>3</sup>	—	—	—
Individual Product Links:	<a href="#">🔗</a>	<a href="#">🔗</a>	<a href="#">🔗</a>	<a href="#">🔗</a>	<a href="#">🔗</a>

— Does Not Apply    🔗 Active Hyperlink  
<sup>1</sup> If OEM model is seamless top, request paper roll holder when ordering  
<sup>2</sup> For OEM models manufactured after 2009 only  
<sup>3</sup> Lift Handles and Top Attachment Knobs for exam table year models 2009 and after are not included

### Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Swatch Card is available from your representative



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2018, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.

Product Hyperlinks Enabled