

PEDIATRIC & YOUTH STRETCHER



THE 360 DEGREE PROTECTION SOLUTION

4160DoST & 462EMPST PEDIATRIC & YOUTH STRETCHERS



INTENDED USE:

The littlest patients need the most protection. Hausted® Pediatric & Youth Stretchers are designed to surround the patients with 360° of safety while maintaining complete access to the patient.

- The Pediatric Stretcher provides safety and security for infants to six years old
- The Youth Stretcher is ideal for patients who are too large for a crib or pediatric stretcher from six years old to adolescence

Pediatric 4160DoST

Cover photo shows optional side and end rail pads and removable, telescoping IV pole (pads: H0001P00B) & (IV: H00018) available at additional cost



4160DoST - Pediatric Product Info

VERSATILITY:

- Dual Side Hydraulic Foot Pedal Controls for easy height adjustment height
- Extra high side rails provide three positions for patient safety and transfer
- Removable end rails to allow easy patient access
- Superior mobility with 8 inch Tente® casters, central brake and steering, and levers conveniently located at all four corners

Youth 462EMPST

Shown with telescoping IV pole (right) (H00018) available at additional cost



462EMPST - Youth Product Info

Specifications:	Pediatric 4160DoST		Youth 462EMPST	
	inch	cm†	inch	cm†
Overall Width:	29.75	75.6	30	76.2
Overall Length:	59.25	150.5	83	210.82
Overall Height: <i>range - minimum</i>	32.25	93.98	28.5	72.39
Overall Height: <i>range - maximum</i>	39	99.06	37	93.98
<i>patient surface</i> Width:	25.25	64.14	25	63.5
<i>patient surface</i> Length:	56	142.24	76	193.04
Spindle Gap: <i>range - minimum</i>	3	7.62	2.5	6.35
Spindle Gap: <i>range - maximum</i>	3.37	8.56	4	10.16
Side Rail Length:	56	142.24	73	185.42
Side Rail Height:	17.75	45.09	16	40.64
End Rail Length:	23	58.42	23	58.42
End Rail Height:	17.75	45.09	16	40.64
<i>easy rolling</i> Tente® Casters:	8	20.3	8	20.3
Fowler:	0 - 45 manual		0 - 90 pneumatic assist	
<i>standard</i> Trendelenburg:	10		12	
<i>reverse</i> Trendelenburg:	10		12	
Weight Capacity: <i>EVENLY DISTRIBUTED</i>	lb	kg†	lb	kg†
	325	147.42	500	226.8

Features:

Hydraulic Controls:	Dual Side Foot Pedals	Dual Side Foot Pedals
---------------------	-----------------------	-----------------------

† standard rounding method used to show metric values as whole numbers

Product Hyperlinks Enabled



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2014, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.