

HORIZON AIRGLIDE™ STRETCHER




YOUR BARIATRIC PATIENT SOLUTION

4E2DPAST / 4D2DPAST HORIZON AIRGLIDE™ BARIATRIC STRETCHER

INTENDED USE:

The Powered Bariatric Stretcher is designed to meet your bariatric patient management needs while versatile enough to be used in every phase of care.



MOBILITY:

- Built-in central steering and braking system with pedals conveniently located at all four corners
- 8 inch Tente® casters allow for smooth movement
- Rechargeable

PATIENT COMFORT:

- Pendant control and caregiver lockout control module for secure patient positioning
- Extra wide 32 inch comfort mattress provides an excellent choice for extended stays
- Full Length, tuck under AirGlide™ Patient Safety Rails with Dual Pneumatic Assists keep patient safe during transport
- Emergency backrest release lever is located conveniently on both sides of unit which allows for quick positioning
- Safe patient egress and reduced risk of back injury for attendant
- Powered unit with head, foot, hi-low and Trendelenburg/Reverse Trendelenburg

Specifications:

	4E2DPAST ¹ 4D2DPAST ¹	
	inch	cm ¹
Overall Width:	35	88.9
Overall Length:	82	208.3
Overall Height: <i>range - minimum</i>	21.5	54.6
Overall Height: <i>range - maximum</i>	29.5	74.9
<i>patient surface</i> Width:	32	81.3
<i>patient surface</i> Length:	77	195.6
<i>AirGlide rail</i> Length:	75	190.5
<i>AirGlide rail</i> Height:	15	38.1
<i>4-way stretch comfort</i> Mattress Pad:	5	12.7
<i>easy rolling</i> Tente® Casters:	8	20.3
	degree of movement ¹	
Fowler: <i>movement</i>	0 - 80	
Knee Gatch: <i>movement</i>	0 - 50	
Trendelenburg: <i>standard</i>	0 - 12	
Trendelenburg: <i>reverse</i>	0 - 12	
Weight Capacity: <i>EVENLY DISTRIBUTED</i>	lb	kg ¹
Shipping Weight: ^{4E2DPAST}	720	327
Shipping Weight: ^{4D2DPAST}	462	209.6
Shipping Weight: ^{4D2DPAST}	440	199.6



Features:

Operation:	Electric/Battery
120/230 volt Option: ¹	<input checked="" type="checkbox"/>
	<input checked="" type="checkbox"/> Standard <input type="checkbox"/> Optional <input type="checkbox"/> Does Not Apply

[4E2/4D2 hyperlink](#)

Color variations may occur due to monitor or printer settings.