

MOBILE SURGI-STRETCHER



ELECTRIC UNITS FOR OPHTHALMOLOGY

5E8 SURGI-STRETCHER SERIES

The Surgi-Stretcher Series is designed as a mobile surgical platform for EYE/ ENT and plastic surgery procedures.

VERSATILITY:

- Articulating head piece with dual operating knobs allow infinite adjustability and precise microscopic movement of the head section
- Designed with a low-profile head piece allowing site access for both superior (over brow) and lateral (side approach) procedures
- Easy-to-use electric/battery powered backrest and leg flex operated controls are available for height adjustment, backrest, leg section, and Trendelenburg/Reverse Trendelenburg



Microscopic articulation of the head piece is achieved by rotating the control knobs to raise/lower the patient's forehead or chin.

STANDARD FEATURES:

- Four wheel brake and steer caster system operable from all four corners
- Easy to clean thermoformed base cover with O2 tank holder and ample storage space
- 3-Stage Telescoping I.V. pole [permanently mounted at foot]
- Stainless steel Retracto® Side Rails with rail caps
- Dual operating control knobs
- Comfort mattress and head pad
- Easy rolling casters
- Low profile articulating head piece

5E8/5D8 ELECTRIC / BATTERY POWERED STANDARD FEATURES:

- Power hand held controls:
 - Height Adjustment
 - Trendelenburg/Reverse Trendelenburg
 - Backrest Adjustment
 - Knee Flex Adjustment
 - Auto Contouring Adjustment
- Nurse Control Panel at foot-end
- Emergency Backrest Release lever to override electric/battery
- Battery back-up allows power functions when not plugged in



Specifications:

	5E8EYEST ¹ 5E8EYXST ¹	
	inch	cm
Overall Width:	29.75	75.6
Overall Width: <i>EYX models only</i>	34.5	87.6
Overall Length:	88.5	224.8
Overall Height: <i>range - minimum</i>	22	55.9
Overall Height: <i>range - maximum</i>	30	76.2
<i>patient surface</i> Width:	25	63.5
<i>patient surface</i> Width: <i>EYX models only</i>	30	76.2
<i>patient surface</i> Length:	81.5	207
<i>retracto® rail</i> Length:	49	124.5
<i>retracto® rail</i> Height:	14	35.6
<i>comfort</i> Mattress Pad:	3	7.6
<i>easy rolling</i> Casters:	8	20.3
	<small>degree of movement[†]</small>	
Backrest: <i>movement</i>	0 - 80	
Knee Flex: <i>movement</i>	0 - 60	
Trendelenburg: <i>standard</i>	10	
Trendelenburg: <i>reverse</i>	10	
Weight Capacity: <small>EVENLY DISTRIBUTED</small>	lb	kg [†]
	720	327
Shipping Weight: <i>EYEST</i>	440 200	
Shipping Weight: <i>EYXST</i>	465 211	

Features:

Operation:	Electric/Battery
<i>dual side hydraulic</i> Foot Pedal Controls: ²	—
<i>foot-end hand crank controlled</i> Backrest:	—
<i>foot-end hand crank controlled</i> Leg Flex:	—
DC/220 volt Option: ¹	■

■ Standard — Not Applicable/Available

¹ same specs for both AC110 volt and the DC220 volt versions

² foot pedal controls adjust height & trendelenburg positioning

[†] standard rounding method used to show metric values as whole numbers

5E8 Product Info



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.