

# MOBILE SURGI-STRETCHER



ELECTRIC UNITS FOR OPHTHALMOLOGY

## 5E82 ELECTRIC SURGI-STRETCHER SERIES

The Surgi-Stretcher Series is designed as a mobile surgical platform for EYE/ENT and plastic surgery procedures.

### VERSATILITY:

- Articulating head piece with dual operating knobs allow infinite adjustability and precise microscopic movement of the head section
- Designed with a low-profile head piece allowing site access for both superior (over brow) and lateral (side approach) procedures
- Battery powered stretcher with easy to use interactive hand pendant.
- Leg flex operated controls are available for height adjustment, backrest, leg section, and Trendelenburg/Reverse Trendelenburg



### STANDARD FEATURES:

- Four wheel brake and steer with electric casters system
- Base cover with built in O2 tank holder and ample storage space
- 3-Stage Telescoping IV pole [permanently mounted at foot]
- Stainless steel Retracto® Side Rails with rail caps
- Low-profile articulating head piece with dual operating control knobs
- Comfort mattress and head pad
- Hand held controls:
  - Height Adjustment
  - Trendelenburg/Reverse Trendelenburg
  - Backrest Adjustment
  - Knee Flex Adjustment
  - Auto Contouring Adjustment
- Emergency Backrest Release to override electric/battery
- Battery back-up allows power functions when not plugged in

### Specifications: 5E82EYEST 5E82EYXST

	standard width		wide width	
	inch	cm†	inch	cm†
Overall Width:	32.00	81.28	37.00	93.98
Overall Length:	92.34	234.54	92.34	234.54
Overall Height: <small>minimum</small>	22.00	55.88	22.00	55.88
Overall Height: <small>maximum</small>	40.00	101.60	40.00	101.60
Patient Surface: <small>width</small>	25.00	63.50	30.00	76.20
Patient Surface: <small>length</small>	85.38	216.87	85.38	216.87
Side Rail Length:	49.00	124.46	49.00	124.46
Side Rail Height:	14.00	35.56	14.00	35.56
Mattress Pad: <small>thickness</small>	3.00	7.62	3.00	7.62
Casters:	5.90	14.99	5.90	14.99
Backrest: <small>movement</small>	degrees [°] 0 - 85		degrees [°] 0 - 85	
Knee Flex: <small>movement</small>	degrees [°] 0 - 60		degrees [°] 0 - 60	
Trendelenburg: <small>standard</small>	15		15	
Trendelenburg: <small>reverse</small>	15		15	
Volts (V)	power 100 - 240		power 100 - 240	
Hertz (Hz)	~ 50/60		~ 50/60	
Weight Capacity: <sup>1</sup> <small>EVENLY DISTRIBUTED</small>	lb	kg†	lb	kg†
	800	362.87	800	362.87
Shipping Weight:	420	190.51	460	208.65

■ Standard — Not Applicable/Available

† standard rounding method used to show metric values as whole numbers



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2023, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.

**graham**  
**field**