

# APC ALL PURPOSE CHAIR



SINGLE PATIENT PLATFORM FOR MULTIPLE PROCEDURES

## STRETCHEROLOGY: APC250ST HYDRAULIC ALL PURPOSE CHAIR

The Surgi-Stretcher Series is designed as a mobile surgical platform for EYE/ENT and plastic surgery procedures.

### INTENDED USE:

The Hausted® All-Purpose Procedure Chair is uniquely designed to provide a single patient platform for a wide variety of procedures in the:

- Emergency/Triage
- Surgery Center/Outpatient Surgery
- Physician's Office
- GI Lab/Endoscopy



**APC250ST**  
HYDRAULIC MODEL  
shown in supine position



[APC Product Info](#)

### VERSATILITY:

- Streamlined handling through treatment, transport and recovery
- Multi-positional, moving from chair to stretcher position in seconds
- Access to steer and brake pedals on all four corners

### MOBILITY:

- Easy to clean thermoformed base with space spacing saving footprint makes chair easy to transport
- Four wheel brake and steer system with 5 inch (12.7 cm) Easy Roll casters

### Specifications: APC250ST

	inch	cm <sup>†</sup>
Overall Width:	27.75	70.5
Overall Length:	75	191
Overall Height: <i>range - minimum</i>	23.5	60
Overall Height: <i>range - maximum</i>	31.5	80
<i>patient surface</i> Width:	22	56
<i>patient surface</i> Length:	70.5	179.1
<i>side rail</i> Length:	29	74
<i>side rail</i> Height:	12	30
<i>comfort</i> Mattress Pad:	3 <sup>1</sup>	7.6 <sup>1</sup>
<i>easy rolling</i> Casters:	5	13
	<small>degree of movement °</small>	
Backrest: <i>movement</i>	0 - 80	
Knee Flex: <i>movement</i>	0 - 75	
Trendelenburg: <i>standard</i>	0 - 10	
Trendelenburg: <i>reverse</i>	0 - 10	
	lb	kg <sup>†</sup>
Weight Capacity: <small>EVENLY DISTRIBUTED</small>	325	147
Shipping Weight:	450	204

**Features:**  
Operation: Hydraulic

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company  
**MDSAP**  
Intertek Medical Device Single Audit Program



+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)

<sup>†</sup> standard rounding method used to show metric values as whole numbers

Product Hyperlinks Enabled