

INTRODUCING THE NEW
HSS SERIES SURGICAL STOOL



HSS500-S SADDLE SEAT STYLE



HSS500-WD WATER DROP SEAT STYLE



HSS500-W WEDGE SEAT STYLE

BUILT FOR PERFORMANCE
AND COMFORT

HAUSTED SURGICAL STOOL

IDEAL FOR EYE, ENT AND PLASTIC PROCEDURES




360° Rotating Urethane Arm Rest for Comfortable Arm Positioning

Adjustable Contoured Backrest

300 POUNDS
300 lb (136 kg)
Weight Capacity
EVENLY DISTRIBUTED

Dual-density Foam Seat with **DURA-GEL®** Insert for Enhanced Comfort



180° Arm Rest Support for more Stable Arm Positioning

Adjustable Height Range
21 - 28 inches (53 - 71 cm)

OPTIONAL ACCESSORY:
Riser Kit Increases Height Range an Additional 2 inches (5 cm)
n No: HSS-Riser

ABS Thermoformed Hood for Easy Cleaning



FOOT-CONTROLLED HYDRAULICS with Three-Spindle Base and Four 3 inch (8 cm) Dual Tente® Casters for Greater Stability



Front Locking System

Complete the HSS Series Surgical Stool Experience with the new Hausted **ESC² Ophthalmology Platform**
Item No: ESC2EYEST

ENGINEERED TO PROVIDE
ENHANCED
SURGEON COMFORT

WARNINGS: Cancer and Reproductive Harm -
www.P66Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



EC REP

CE Partner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.eu



© 2020, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE Casters, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com

