

TREATMENT TABLES





Special Features

Our pass-through drawers can be opened from either side of the table, so you'll always have easy access to supplies no matter where your table sits.



Treatment Table, Flat Top

4100-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:
18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 170 lbs / 77 kg

Treatment Table with Adjustable Back and Contour Top

4105-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:
18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 190 lbs / 86 kg

Treatment Table, Flat Top, Shelf and Two Drawers

4103-DD-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:
18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 170 lbs / 77 kg

Treatment Table with Adjustable Back, Contour Top, Shelf and Two Drawers

4105-DD-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:
18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 210 lbs / 95 kg

Standard Color Options(-xx)

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Note: Due to the many variations in monitors and browsers, color samples may appear different. Swatch Card is available from your representative.

1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can
be found on our website.
www.grahamfield.com

Follow us on



© 2018, GF Health Products, Inc. All Rights Reserved.
Hausted is a trademark of GF Health Products, Inc. GF Health
Products, Inc. is an ISO 13485 Certified Company.
GF1800081RevC19

graham
field