



SUPPORT SURFACES




# HAUSTED PROCEDURE CHAIRS REPLACEMENT PADS & FABRICS

The Hausted® Clinical Procedure Chair Pads are designed to provide excellent pressure reduction, comfort, support, and durability.



## STANDARD PROCEDURE CHAIR STRETCH FABRIC








049 BLUE

### PRIMASTRETCH

Primastretch is engineered with high tensile strength for durability. The 4-way stretch, urethane coated, nylon cover is breathable, reduces friction, and skin shear while helping to keep the patient cool and comfortable.

## REPLACEMENT PAD SETS:


Procedure Chair Model N°:	049 BLUE	051 NAVY	050 BLACK	037 BLACK	043 SPRUCE
	PRIMASTRETCH FABRIC			SURE-CHEK FUSION III	
ESC2EYEST	 H080901-049	 H080901-051	 H080901-050	 H080901-037	 H080901-043
APC250ST	HP150830384-049	HP150830384-051	HP150830384-050	HP150830384-037	HP150830384-043
EPC500ST	H081185-049	H081185-051	H081185-050	H081185-037	H081185-043
MBC000ST	H134604-049	H134604-051	H134604-050	H134604-037	H134604-043
VIC429ST	H134602-049	H134602-051	H134602-050	H134602-037	H134602-043
<b>Fabric Color:</b>	<b>Standard</b>	<b>Optional Fabric</b>		<b>Upgraded Fabric</b>	

swatches are digital reproductions, perceived texture and coloring may vary on finished product, color variations may occur due to monitor or printer settings.



**REPLACEMENT PAD SETS: Optional Fabric**

Procedure Chair Model N°:	427 BLUE-RIDGE	8602 OATMEAL	8604 CONCRETE	8567 DOVE	454 ROYAL BLUE	432 IMPERIAL BLUE	8584 TAUPE	8516 CHESTNUT	8615 BERRY
	<b>SPRADLING PATRIOT PLUS VINYL</b> <span style="border: 1px solid black; padding: 2px;">NEO</span>								
ESC2EYEST	H080901-427	H080901-8602	H080901-8604	H080901-8567	H080901-454	H080901-432	H080901-8584	H080901-8516	H080901-8615
APC250ST	HP150830384-427	HP150830384-8602	HP150830384-8604	HP150830384-8567	HP150830384-454	HP150830384-432	HP150830384-8584	HP150830384-8516	HP150830384-8615
EPC500ST	H081185-427	H081185-8602	H081185-8604	H081185-8567	H081185-454	H081185-432	H081185-8584	H081185-8516	H080901-8516
MBC000ST	H134604-427	H134604-8602	H134604-8604	H134604-8567	H134604-454	H134604-432	H134604-8584	H134604-8516	H134604-8615
VIC429ST	H134602-427	H134602-8602	H134602-8604	H134602-8567	H134602-454	H134602-432	H134602-8584	H134602-8516	H134602-8615



037 BLACK



043 SPRUCE

**STRETCH URETHANE SURE-CHEK® FUSION III**

The Sure-Chek Fusion III is engineered for maximum performance. The 4-way stretch fluid-resistant urethane cover reduces friction and skin shear while helping to keep the patient cool and comfortable. The fabric provides maximum resistance against spills, staining, and odors. The fabric incorporates antimicrobial additives designed to help protect the fabric against odor-forming bacterial growth to reduce undesirable organic odors.




051 NAVY




050 BLACK

**PRIMASTRETCH**


Primastretch is engineered with high tensile strength for durability. The 4-way stretch urethane coated nylon cover is breathable, reduces friction, and skin shear while helping to keep the patient cool and comfortable.



8602 OATMEAL




8604 CONCRETE



8567 DOVE

**PATRIOT PLUS POLYVINYL** NEO


The Patriot-Plus Neo is a polyvinyl material designed with their antimicrobial protection and resistance to abrasion and stains from the demands of the healthcare environment. The fabric is easy to clean and contains anti-static properties.




427 BLUERIDGE




8584 TAUPE




8516 CHESTNUT



432 IMPERIAL BLUE



454 ROYAL BLUE




8615 BERRY

swatches are digital reproductions. perceived texture and coloring may vary on finished product. color variations may occur due to monitor or printer settings.



**INDIVIDUAL REPLACEMENT PADS:**

To order an individual replacement pad, use the item number below replacing xxx with the standard, optional, or upgraded fabric code found on the preceding pages



**APC**

	FOOT		LEG		SEAT		BACK		HEAD	
Overall	HP150830389-xxx		HP150830388-xxx		HP150830387-xxx		HP150830386-xxx		HP150830385-xxx	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Width:	14.00	35.56	22.00	55.88	22.60	57.40	22.00	55.88	20.50	52.07
Length:	7.75	19.68	10.00	25.40	20.00	50.80	20.00	50.80	7.50	19.05
Height:	3.00	7.62	3.00	7.62	3.00	7.62	3.00	7.62	2.75	6.98

**MBC**

	FOOT		LEG		SEAT		BACK		HEAD	
Overall	HP150830389-xxx		HP150830388-xxx		HP150830387-xxx		HP150830386-xxx		HP150830385-xxx	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Width:	14.00	35.56	22.00	55.88	22.60	57.40	22.00	55.88	22.00	55.88
Length:	7.25	18.42	10.00	25.40	20.00	50.80	20.00	50.80	12.50	31.75
Height:	3.00	7.62	3.00	7.62	3.00	7.62	3.00	7.62	3.00	7.62

**ESC<sup>2</sup>**

	FOOT		LEG		SEAT		BACK		HEAD	
Overall	H080862-xxx		H080863-xxx		H080864-xxx		H080865-xxx		H080866-xxx	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Width:	20.50	52.07	25.00	63.50	26.00	66.04	25.00	63.50	10.00	25.40
Length:	7.50	19.05	10.30	26.16	18.90	48.01	21.50	54.61	10.00	25.40
Height:	2.75	6.98	2.75	6.98	2.75	6.98	2.75	6.98	2.25	5.72

**EPC**

	FOOT		LEG		SEAT		BACK	
Overall	H080862-xxx		H080863-xxx		H080864-xxx		H081166-xxx	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Width:	20.50	52.07	25.00	63.50	26.00	66.04	25.00	63.50
Length:	7.50	19.05	10.30	26.16	18.90	48.01	31.25	79.38
Height:	2.75	6.98	2.75	6.98	2.75	6.98	2.75	6.98

**VIC**

	FOOT		LEG		SEAT		BACK-L		BACK-R	
Overall	HP150830389-xxx		HP150830388-xxx		HP150830387-xxx		H13460F-xxx		H13460W-xxx	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†
Width:	14.00	35.56	22.00	55.88	22.60	57.40	13.00	193.04	8.25	20.96
Length:	7.25	18.42	10.00	25.40	20.00	50.80	31.50	80.01	31.50	80.01
Height:	3.00	7.62	3.00	7.62	3.00	7.62	2.00	5.08	2.00	5.08

## ACCESSORY PADS:

Item	Description	Size <sup>4</sup>		UOM <sup>3</sup>	MBC000ST	VIC429ST	APC250ST	ESC2EYEST	EPC500ST
		inch	cm <sup>1</sup>						
<b>Accessory Mounting Rails</b>									
HSA080013-xxx	Shoulder Rail/Extensions <sup>2</sup>	8 x 8	20 x 20	pair				■	■
<b>Armboards</b>									
H065990	Seated without Pad			each	■	■	■	■	■
H0659900-xxx	Seated Assembly with Pad	26 x 6	66 x 15.2	each	■	■	■	■	■
HSA041300-xxx	with Pad and Clamp	2	5.08					■	■
H06123900-xxx	with Pad for Reclined						■		
HSA080029-xxx	Ortho Hand Surgery Board							■	■
<b>Foot Extension</b>									
H080015-xxx	Chair <sup>2</sup>				■	■	■		
<b>Contoured Headrests</b>									
H0103ESC2-xxx	Set of 3 Sizes below <sup>2</sup>			set of 3				■	
H080988-xxx	Individual - small <sup>2</sup>	4	10	each				■	
H080989-xxx	Individual - medium <sup>2</sup>	4.5	11	each				■	
H080990-xxx	Individual - large <sup>2</sup>	5	13	each				■	
<b>Hand Boards</b>									
HSA080029-xxx	Hand Surgery Board <i>-Orthopedic</i>			each				■	■



### Foot Extension

Provides extra leg length 14 inch (35.6cm) to foot end

Item No: H080015-xxx<sup>2</sup>



### Shoulder Rail Extensions

Provides Patient Support and Comfort

Item No: HSA080013-xxx<sup>2</sup>

### Contoured Headrest Kit (3pc):

Includes 3 headrests<sup>2</sup> (small, medium, and large)

Item No: H0101ST-xxx<sup>2</sup>



### Armboard

Mounts to surgical accessory rail for greater flexibility in patient positioning

Item No: H07633400

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved. Hausted is a registered trademark of GF Health Products, Inc.