

## SUPPORT SURFACES



# HAUSTED UPGRADE AND REPLACEMENT STRETCHER PADS & FABRICS

### 600C COMFORT CARE SERIES UPGRADE STRETCHER PADS

The Hausted<sup>®</sup> Clinical Comfort Upgrade Stretcher Pads are the best choice when selecting a non-powered support surface with excellent pressure reduction properties compared to a standard pad.

#### Features:

##### 5 inch thick foam pad

- Tri-laminate foam pad design for better comfort and support
- Contour cut grid allows for better patient immersion, and greater airflow
- Better pressure redistribution properties than standard pad offerings
- Hook attachment on base
- Easy to clean and maintain

#### Cover Features:

- 4-way stretch cover material is durable and provides a softer feel.
- Infection Control is aided by Antimicrobial properties; Bacteria resistant, anti-fungal, and non-allergenic
- Flame Resistant
- Fluid Resistant
- Anti-static
- Stain Resistant



### H0 STANDARD SERIES STRETCHER PADS

The Hausted<sup>®</sup> Standard Stretcher Pads offer exceptional cushioning, support, and pressure redistribution. The fabric is durable and offers Low friction for easy patient transfers.

#### Features:

##### Standard foam pad height

- 2 layer foam pad for comfort and durability
- Hook attachment on base
- Easy to clean and maintain

#### Cover Features:

- Herculite<sup>®</sup> Sure-Chek<sup>™</sup> cover material is comfortable and durable
- Infection Control is aided by Antimicrobial properties; Bacteria resistant, anti-fungal, and non-allergenic
- Flame Resistant
- Fluid Resistant
- Anti-static
- Stain Resistant



**STANDARD SERIES**  
2-layer foam construction for comfort

## STANDARD STRETCHER PAD

The Hausted Stretcher Pad comes standard on our stretchers. The 2-layer foam construction provides patient comfort. The Guardian EcoSense Vinyl fabric cover is fluid resistant, durable and easy to clean.



**600C — COMFORT CARE SERIES**  
3-layer construction with contour-cut grid pattern provides greater airflow and improves patient immersion

## COMFORT CARE STRETCHER UPGRADE PAD

The Hausted Comfort Care Stretcher Pad is an excellent upgrade option when longer patient stays may occur. The 5-inch (13 cm) thickness provides maximum support and pressure reduction properties compared to a standard pad. The contour-cut grid and 3-layer foam design provides greater airflow, improves patient immersion, and extended patient use. The fluid-resistant 4-way stretch urethane (Sure-Chek Fusion III) cover keeps the patient cool and comfortable while providing maximum resistance against spills, staining and odors.

Stretcher Model N°:	Standard Pads			Upgraded Pads
	SENTRY Guardian EcoSense Fabric 1324	BLACK Herculite® Lectrolite™ Fabric 040	BLACK Sure-Chek® Comfort Fabric 036	SPRUCE Sure-Chek® Fusion III Fabric <sup>2</sup> 043 047 036 NAVY Sure-Chek® Fusion III Fabric <sup>2</sup> 036 BLACK Sure-Chek® Fusion III Fabric <sup>2</sup>
4160D0ST	■ H031612B	■ H031612-040	■ H031612-036	—
462EMPST	■ H066283B	■ H06189360-040	■ H06189360-036	■ 600C-462EMPST-xxx
462EMCST	■ H06189360B	■ H066283-040	■ H066283-036	■ 600C-462EMCST-xxx
462FMCST	■ H06189460B	■ H06189460-040	■ H06189460-036	■ 600C-462FMCST-
462GMCST	■ H06189360B	■ H06189360-040	■ H06189360-036	■ 600C-462GMCST-xxx
462HMCST	■ H06189460B	■ H06189460-040	■ H06189460-036	■ 600C-462HMCST-xxx
462APAST	■ H06717860B	■ H06717860-040	■ H06717860-036	■ 600C-462APAST-xxx
462CPAST	■ H06717860B	■ H06717860-040	■ H06717860-036	■ 600C-462EMCST-xxx
462DPAST	■ H06896560B	■ H06896560-040	■ H06896560-036	■ 600C-462DPAST-xxx
493RPCST	■ H03165560B	■ H03165560-040	■ H03165560-036	—
578EYEST	■ H075112B ★H069538B	■ H075112-040 ★H069538-040	■ H075112-036 ★H069538-036	■ 600C-578EYEST-xxx
578EYXST	■ H075255B ★H069538B	■ H075255-040 ★H069538-040	■ H075255-036 ★H069538-036	■ 600C-578EYXST-xxxx
5E82EYEST	■ H081235B	■ H081235-040	■ H081235-036	■ 600C-5E82EYEST-xxx
5E82EYXST	■ H081328B	■ H081328-040	■ H081328-036	■ 600C-5E82EYXST-xxx


Bariatric Stretcher Model N°:	Upgraded		
	NAVY Primastretch Fabric 051	SPRUCE Sure-Chek® Fusion III Fabric <sup>2</sup> 043 047 036	NAVY Sure-Chek® Fusion III Fabric <sup>2</sup> 036
5B800ST	■ H081385-051	■ 600C-081385-xxx	—
5B800ST-PD	■ H081385-051	■ 600C-081385-xxx	—

**Model Number Shown:**  
orders frame with standard pad and standard fabric

**Ordering Options:**  
add an 'x' at the end of the model number (orders frame only), then order desired pad and fabric as a separate item


- Standard Fabric Color
- Optional Fabric Color
- Upgraded Fabric & Pad
- ★ Matching Head Pad Set
- Not Available

<sup>2</sup> multiple fabrics options available  
- replace 'xxx' with color code



## STRETCHER FABRICS


The **Hausted** stretcher pads use fabrics designed to meet the demands in a healthcare environment. Pads are fluid-resistant and resist staining from most contaminants such as medications, oils, and more. The fabric contains antimicrobial additives designed to help protect the fabric against odor-forming bacterial growth and reduces undesirable organic odors.



1324 **SENTRY**

### STANDARD STRETCHER PAD FABRICS GUARDIAN ECOSENSE FABRIC

**Guardian EcoSense** is a polyvinyl material designed with antimicrobial protection and resistance to abrasion and stains. The fabric is easy to clean and contains anti-static properties.




036 **BLACK**

### STANDARD STRETCHER PAD OPTIONAL FABRICS

#### SURE-CHEK® COMFORT FABRIC

**Sure-Chek Comfort** is a fluid and moisture-resistant cover known for its durability and easy cleaning. It is a vinyl laminate material that conforms to the body and allows the inner construction of the stretcher pad to provide comfort and support.



040 **BLACK**

#### HERCULITE® LECTROLITE COMFORT VINYL

**Lectrolite** is a conductive vinyl laminated fabric to help control electrostatic discharge. This conductive fabric cover also prevents the accumulation of a static charge that can interfere with sensitive electronic equipment. Lectrolite conforms comfortably to the body thereby allowing the patient to benefit from the comfort features of the stretcher pad.



051 **NAVY**

### BARIATRIC STANDARD STRETCHER PAD FABRIC PRIMASTRETCH FABRIC

**Primastretch** is engineered with high tensile strength for durability for bariatric use. The 4-way stretch Urethane coated nylon cover is breathable and reduces friction and skin shear while helping to keep the patient cool and comfortable.



043  
**SPRUCE**



047  
**NAVY**



036  
**BLACK**

### UPGRADED STRETCHER PAD FABRICS

#### SURE-CHEK® FUSION III®

The **Sure-Chek Fusion III** is engineered for maximum performance. The 4-way stretch, fluid-resistant, urethane cover reduces friction and skin shear while helping to keep the patient cool and comfortable. The fabric provides maximum resistance against spills, staining, and odors. The fabric incorporates antimicrobial additives designed to help protect the fabric against odor-forming bacterial growth to reduce undesirable organic odors.

swatches are digital reproductions, perceived texture and coloring may vary on finished product, color variations may occur due to monitor or printer settings.

## HAUSTED STRETCHER REPLACEMENT PADS FOR ACCESSORIES

Item	Description	UOM	Width		Length		Height/Thickness		4160DOST	462EMPST	462EMCST	462FMCST	462GMCST	462HMCST	462APAST	462CPAST	462DPAST	493RPCST	5B80OST	5B80OST-PD	5E8EYEST	5D8EYEST	5E8EYXST	5D8EYXST	578EYEST	578EYXST	5E82EYEST	5E82EYXST	
			inch	cm <sup>†</sup>	inch	cm <sup>†</sup>	inch	cm <sup>†</sup>																					
<b>Armboard Pads</b>																													
HP150830167-xxx	Surgical	1 ea	6.00	15.24	26.00	66.04	2.25	5.72																					
<b>Handboard Pads</b>																													
H081315-xxx	Ortho Hand Surgery Board	1 ea	12.50	31.75	27.10	68.83	2.25	5.72																					
<b>Headrests</b>																													
H0101ST-xxx	Set of 3 of below. <sup>2</sup>	1 set																											
HSA063500-xxx	- large <sup>2</sup>	1 ea	10.50	26.67	9.50	24.13	4.00	10.16																					
HP150830448-xxx	- medium <sup>2</sup>	1 ea	10.50	26.67	9.50	24.13	3.00	7.62																					
HP150830447-xxx	- small <sup>2</sup>	1 ea	10.50	26.67	9.50	24.13	2.00	5.08																					
<b>Rail Cover Pads</b>																													
HPD0610-xxx	Retracto® Pad 4 Leg Set	1 pair	1.25	3.18	48.50	123.19	15.00	38.10																					
HPD0612-xxx	Retracto® Pad 6 Leg Set	1 pair	1.25	3.18	54.00	137.16	12.00	30.48																					
H00011P00-xxx	Side & End Pad Set of below.	set of 4 (2 ea)																											
H030158-xxx	Side Rail Pad	1 ea	1.00	2.54	55.12	140.00	15.50	39.37																					
H030159-xxx	End Rail Pad	1 ea	1.00	2.54	23.00	58.42	15.50	39.37																					
H00P11R00-xxx	Side & End Pad Set of below.	set of 4 (2 ea)																											
H066045-xxx	Side Rail Pad	1 ea	1.00	2.54	73.00	185.42	12.25	31.12																					
H030159-xxx	End Rail Pad	1 ea	1.00	2.54	23.00	58.42	13.00	33.02																					

<sup>†</sup> standard rounding method used    <sup>2</sup> multiple fabrics options available - replace 'xxx' with color code - see following pages    <sup>8</sup> approximate measurement



### Rail Covers

Slip-on side and end covers for extra patient protection

Item No: H00011P00-xxx<sup>2</sup>

### Contoured Headrest Kit (3pc):

Includes 3 headrests<sup>2</sup> (small, medium, and large)

Item No: H0101ST-xxx<sup>2</sup>



### Armboard

Mounts to surgical accessory rail for greater flexibility in patient positioning

Item No: H07633400

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2016, GF Health Products, Inc. All Rights Reserved. Hausted is a registered trademark of GF Health Products, Inc.