

ELECTRIC ALL PURPOSE
BARIATRIC STRETCHER



5B800ST ELECTRIC ALL PURPOSE BARIATRIC STRETCHER

The powered all purpose bariatric stretcher is designed to meet your patient management needs while versatile enough to be used in every phase of care.



- **Optional Navigator SmartDrive System:** fully powered drive system for stretcher positioning and patient transport - forward, backward, and side to side **5B800ST-PD**
- **800 lb (362.87 kg) Weight Capacity** – EVENLY DISTRIBUTED¹
- Extra-wide 30 inch (76 cm) patient surface provides greater patient comfort for extended stays
- Stainless steel Retracto[®] Side Rails with aluminum power coated rail caps
- Integrated Push Handles for Increased Staff Safety
- Electric caster system allows for easy maneuverability and transfer of patients
- Lithium battery allows for increased battery life
- Comfort Mattress

♦ See Chart for Specifications

5B800ST Product Info



STANDARD FEATURES:

- Four wheel Brake and Steer with Electric caster system allows for easy maneuverability
- Easy to clean thermoformed base cover with built in O2 tank holder and ample storage space
- 3-Stage Telescoping IV pole [permanently mounted at head]
- Pendant controls:
 - Height Adjustment
 - Trendelenburg/Reverse Trendelenburg
 - Backrest Adjustment
 - Knee Flex Adjustment
 - Auto Contouring Adjustment
- Emergency Backrest Release to override electric/battery
- Battery back-up allows power functions when not plugged in



Model: 5B800ST
bariatric

	inch	cm [†]
Overall Width:	37.00	93.98
Overall Length:	88.38	224.60
Overall Height: <i>minimum</i>	22.00	55.88
Overall Height: <i>maximum</i>	40.00	101.60
Patient Surface: <i>width</i>	30.00	76.20
Patient Surface: <i>length</i>	76.44	194.16
Side Rail Length:	62.60	159.00
Side Rail Height:	14.00	35.56
Mattress Pad: <i>thickness</i>	5.00	12.70
Caster Diameter:	5.90	14.99

	features	
Rail Type:	Retracto™	
Caster Type:	dual surface	
	degrees (+)	
Backrest Movement:	0 - 85	
Knee Flex Movement:	0 - 60	
Trendelenburg: <i>standard</i>	15	
Trendelenburg: <i>reverse</i>	15	
	power	
Volts (V)	100 - 240	
Hertz (Hz)	~ 50/60	
	lb	kg [†]
Weight Capacity: ¹ EVENLY DISTRIBUTED	800	362.87
Shipping Weight: ⁸	500	226.80

¹ maximum weight capacity - evenly distributed - refer to user manual
⁸ approximate measurement † standard rounding method used



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2024, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.